


# ANTIMICROBIAL REUSABLE AFFORDABLE FILTER FOR FACE MASKS

**BreaSAFE®**  
**FILTER CLASSIC**  
10 × 15

NANOFIBROUS FILTER PAD  
FOR POCKET FACE MASKS



1 package contains 10 pcs (10 × 15 cm)

\*Other dimensions at the customer's request for B2B purchases

## **BreaSAFE® FILTER transforms basic face masks into an effective weapon against viruses and bacteria**

Unique Czech nanofibrous membrane thanks to active silver effectively eliminates microorganisms, and therefore significantly increases protective factor of commonly available pocket face masks (purchased or home-made).

Super-thin fibers and small pores ensure high effectiveness at capturing particles, bacteria, viruses and allergens.

### **Instructions for use:**


1. Wash your hands thoroughly.
2. Remove BreaSAFE® FILTER from its package and insert it into the pocket of the face mask
3. Place the face mask on your face following instructions of the manufacturer.

### **Instructions for maintenance:**

BreaSAFE® FILTER can be:

- hand washed in disinfectant based on alcohol or 2-propanol
- sterilized in an autoclave or an oven (max. 70 °C)
- sterilized by ironing (at minimum temperature)

## 97% EFFICIENCY



*Proudly developed and made  
in the Czech Republic.*

[www.BreaSAFE.com](http://www.BreaSAFE.com) | [www.nano4fibers.com](http://www.nano4fibers.com)



**NANO4FIBERS**  
**PARDAM**

SIMPLE BUT EFFECTIVE PROTECTION