


ANTIMICROBIAL REUSABLE AFFORDABLE FILTER FOR FACE MASKS

BreaSAFE®
FILTER MAX
15 × 20

NANOFIBROUS FILTER PAD
FOR POCKET FACE MASKS



1 package contains 10 pcs (15 × 20 cm)

*Other dimensions at the customer's request
for B2B purchases

BreaSAFE® FILTER transforms basic face masks into an effective weapon against viruses and bacteria

Unique Czech nanofibrous membrane thanks to active silver effectively eliminates microorganisms, and therefore significantly increases protective factor of commonly available pocket face masks (purchased or home-made).

Super-thin fibers and small pores ensure high effectiveness at capturing particles, bacteria, viruses and allergens.

Instructions for use:


1. Wash your hands thoroughly.
2. Remove BreaSAFE® FILTER from its package and insert it into the pocket of the face mask
3. Place the face mask on your face following instructions of the manufacturer.

Instructions for maintenance:

BreaSAFE® FILTER can be:

- hand washed in disinfectant based on alcohol or 2-propanol
- sterilized in an autoclave or an oven (max. 70 °C)
- sterilized by ironing (at minimum temperature)

97% EFFICIENCY



*Proudly developed and made
in the Czech Republic.*

www.BreaSAFE.com | www.nano4fibers.com



NANO4FIBERS
PARDAM

SIMPLE BUT EFFECTIVE PROTECTION